

XenoWorks Analog Injecor (BRI) Syringe Volume Chart

<b>Syringe Part Number</b>	<b>Syringe Volume</b>	<b>Volume Per Revolution Course Dial</b>	<b>Volume Per Revolution Fine Dial</b>
V001180	10µl	162nl	16.2nl
V001181	25µl	325nl	32.5nl
V001182	50µl	650nl	65.0nl
V001183	100µl	1.3µl	130nl
V001184	250µl	2.6µl	0.26µl
V001185	500µl	5.2µl	0.52µl
V001186	1000µl	10.4µl	1.04µl