



Octagonal Chamber Platform for Warner Series 20 (Warner P and PH)



Used with following Sutter Instrument Stages

MP-78	MPC-78	MT-78-FS
MT-1078	MT-2078	MT-2278

Click on product above to view detailed information

Part No	Material	Overall	ID	OD
X040505	Aluminum	101.6x152.4mm(4x6in)	53.3mm(2.1in)	58.4mm(2.3in)
X040507	Stainless St	101.6x152.4mm(4x6in)	53.3mm(2.1in)	58.4mm(2.3in)

